

## KIRSCHNER WIRE DRILLS

trocar-trocar; 2.5 / 120mm

box with 10 pieces



**SKU:** 212-524-25

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>SPECIFICATIONS</b>	<a href="#">trocar-trocar</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>SPECIFICATIONS</b>	<a href="#">trocar-trocar</a>