

KIRSCHNER WIRE DRILLS

trocar-trocar; 1.0 / 160mm

box with 10 pieces



SKU: 212-526-10

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.10.80
SPECIFICATIONS	trocar-trocar

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.10.80
SPECIFICATIONS	trocar-trocar