

KIRSCHNER WIRE DRILLS

trocar-trocarr; 2.2 / 160mm

box with 10 pieces



SKU: 212-526-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.10.92
SPECIFICATIONS	trocar-trocar

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
Medicon	76.10.92
SPECIFICATIONS	trocar-trocar