

KIRSCHNER WIRE DRILLS

trocar-trocar; 1.6 / 310mm

box with 10 pieces



SKU: 212-528-16

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX957
AUTHORS	KIRSCHNER
Medicon	76.13.16
SPECIFICATIONS	trocar-trocar

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX957
AUTHORS	KIRSCHNER
Medicon	76.13.16

SPECIFICATIONS

[trocar-trocar](#)