

KIRSCHNER WIRE DRILLS

troc-ar-round; 1.4 / 140mm

box with 10 pieces



SKU: 212-505-14

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
SPECIFICATIONS	troc-ar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
AUTHORS	KIRSCHNER
SPECIFICATIONS	troc-ar-round