

KIRSCHNER WIRE DRILLS

trocar-round; 1.5 / 160mm

box with 10 pieces



SKU: 212-506-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX135
KLS Martin	22-631-15
AUTHORS	KIRSCHNER
Medicon	76.10.65
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX135
KLS Martin	22-631-15

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.65

SPECIFICATIONS

[trocar-round](#)