

## KIRSCHNER WIRE DRILLS

trocar-round; 2.0 / 160mm

box with 10 pieces



**SKU:** 212-506-20

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX140
<b>KLS Martin</b>	22-631-20
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.10.70
<b>SPECIFICATIONS</b>	<a href="#">trocar-round</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX140
<b>KLS Martin</b>	22-631-20

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.10.70

**SPECIFICATIONS**

[trocar-round](#)