

KIRSCHNER WIRE DRILLS

trocar-round; 2.2 / 160mm

box with 10 pieces



SKU: 212-506-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX142
KLS Martin	22-631-22
AUTHORS	KIRSCHNER
Medicon	76.10.72
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX142
KLS Martin	22-631-22

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.72

SPECIFICATIONS

[trocar-round](#)