

## KIRSCHNER WIRE DRILLS

trocar-round; 1.5 / 310mm

box with 10 pieces



**SKU:** 212-508-15

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX035
<b>KLS Martin</b>	22-633-15
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.12.65
<b>SPECIFICATIONS</b>	<a href="#">trocar-round</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX035
<b>KLS Martin</b>	22-633-15

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.12.65

**SPECIFICATIONS**

[trocar-round](#)