

KIRSCHNER WIRE DRILLS

trocar-round; 2.0 / 310mm

box with 10 pieces



SKU: 212-508-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX040
KLS Martin	22-633-20
AUTHORS	KIRSCHNER
Medicon	76.12.70
SPECIFICATIONS	trocar-round

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX040
KLS Martin	22-633-20

AUTHORS

[KIRSCHNER](#)

Medicon

76.12.70

SPECIFICATIONS

[trocar-round](#)