

KIRSCHNER WIRE DRILLS

trocar-round; 2.2 / 310mm

box with 10 pieces



SKU: 212-508-22

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX042 |
| KLS Martin | 22-633-22 |
| AUTHORS | KIRSCHNER |
| Medicon | 76.12.72 |
| SPECIFICATIONS | trocar-round |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX042 |
| KLS Martin | 22-633-22 |

AUTHORS

[KIRSCHNER](#)

Medicon

76.12.72

SPECIFICATIONS

[trocar-round](#)