

KIRSCHNER WIRE DRILLS

trocar-flat; 1.8 / 140mm

box with 10 pieces



SKU: 212-515-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| SPECIFICATIONS | trocar-flat |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| SPECIFICATIONS | trocar-flat |