

## KIRSCHNER WIRE DRILLS

trocar-flat; 1.6 / 160mm

box with 10 pieces



**SKU:** 212-516-16

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX116
<b>KLS Martin</b>	22-630-16
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.10.16
<b>SPECIFICATIONS</b>	<a href="#">trocar-flat</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX116
<b>KLS Martin</b>	22-630-16

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.10.16

**SPECIFICATIONS**

[trocar-flat](#)