

KIRSCHNER WIRE DRILLS

trocar-flat; 1.8 / 160mm

box with 10 pieces



SKU: 212-516-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX118 |
| KLS Martin | 22-630-18 |
| AUTHORS | KIRSCHNER |
| Medicon | 76.10.18 |
| SPECIFICATIONS | trocar-flat |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| Aesculap | LX118 |
| KLS Martin | 22-630-18 |

AUTHORS

[KIRSCHNER](#)

Medicon

76.10.18

SPECIFICATIONS

[trocar-flat](#)