

## KIRSCHNER WIRE DRILLS

trocar-flat; 2.5 / 160mm

box with 10 pieces



**SKU:** 212-516-25

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX125
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.10.25
<b>SPECIFICATIONS</b>	<a href="#">trocar-flat</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX125
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.10.25

## SPECIFICATIONS

[trocar-flat](#)