

## KIRSCHNER WIRE DRILLS

trocar-flat; 1.2 / 310mm

box with 10 pieces



**SKU:** 212-518-12

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">Box</a>             |
| <b>TYPE</b>           | <a href="#">Kirschner wires</a> |
| <b>Aesculap</b>       | LX012                           |
| <b>KLS Martin</b>     | 22-632-12                       |
| <b>AUTHORS</b>        | <a href="#">KIRSCHNER</a>       |
| <b>Medicon</b>        | 76.15.12                        |
| <b>SPECIFICATIONS</b> | <a href="#">trocar-flat</a>     |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                                 |
|-------------------|---------------------------------|
| <b>UNIT</b>       | <a href="#">Box</a>             |
| <b>TYPE</b>       | <a href="#">Kirschner wires</a> |
| <b>Aesculap</b>   | LX012                           |
| <b>KLS Martin</b> | 22-632-12                       |

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.15.12

**SPECIFICATIONS**

[trocar-flat](#)