

KIRSCHNER WIRE DRILLS

trocar-flat; 2.0 / 310mm

box with 10 pieces



SKU: 212-518-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX020
KLS Martin	22-632-20
AUTHORS	KIRSCHNER
Medicon	76.15.20
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX020
KLS Martin	22-632-20

AUTHORS

[KIRSCHNER](#)

Medicon

76.15.20

SPECIFICATIONS

[trocar-flat](#)