

## KIRSCHNER WIRE DRILLS

trocar-flat; 2.2 / 310mm

box with 10 pieces



**SKU:** 212-518-22

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX022
<b>KLS Martin</b>	22-632-22
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	76.15.22
<b>SPECIFICATIONS</b>	<a href="#">trocar-flat</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">Box</a>
<b>TYPE</b>	<a href="#">Kirschner wires</a>
<b>Aesculap</b>	LX022
<b>KLS Martin</b>	22-632-22

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

76.15.22

**SPECIFICATIONS**

[trocar-flat](#)