

KIRSCHNER WIRE DRILLS

trocar-flat; 2.5 / 310mm

box with 10 pieces



SKU: 212-518-25

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX025
KLS Martin	22-632-25
AUTHORS	KIRSCHNER
Medicon	76.15.25
SPECIFICATIONS	trocar-flat

ZUSÄTZLICHE INFORMATIONEN

UNIT	Box
TYPE	Kirschner wires
Aesculap	LX025
KLS Martin	22-632-25

AUTHORS

[KIRSCHNER](#)

Medicon

76.15.25

SPECIFICATIONS

[trocar-flat](#)