

KIRSCHNER WIRE DRILLS

trocar-trocar; 1.5 / 60mm

box with 10 pieces



SKU: 212-521-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| SPECIFICATIONS | trocar-trocar |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | Box |
| TYPE | Kirschner wires |
| AUTHORS | KIRSCHNER |
| SPECIFICATIONS | trocar-trocar |