

KLEINERT-KUTZ TENDON HOLDING FORCEPS

2.0mm; straight

20.5cm (8.2")




SKU: 195-248-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-681-20
AUTHORS	KLEINERT-KUTZ
Medicon	78.14.40
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-681-20
AUTHORS	KLEINERT-KUTZ
Medicon	78.14.40

SPECIFICATIONS

[straight](#)