

KLEPPINGER ELECTRODE; BIPOAR; INSERT

serrated; for Ø 5.0mm shaft

for 628-201-53 and 628-205-53



SKU: 628-248-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes