

## KNOT PUSHER; FENESTRATED

axial handle curved; with ratchet

Luer-Lock; Ø 5.0mm; 330mm



**SKU:** 612-295-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Suture Sets](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Suture Sets](#)