


KOCHER-DELICATE HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

16.0cm (6.4")



SKU: 040-352-16

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH176 |
| KLS Martin | 13-384-16 |
| AUTHORS | DELICATE , KOCHER |
| Medicon | 16.12.86 |
| SPECIFICATIONS | serrated , straight , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH176 |
| KLS Martin | 13-384-16 |

AUTHORS

[DELICATE](#), [KOCHER](#)

Medicon

16.12.86

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)