

KOCHER INTESTINAL CLAMP

curved; serrated longitudinal; elastic

28.0cm (11.2")

SKU: 160-235-28

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA157
KLS Martin	30-305-28
AUTHORS	KOCHER
Medicon	50.23.28
SPECIFICATIONS	curved , serrated longitudinal



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA157
KLS Martin	30-305-28

AUTHORS

[KOCHER](#)

Medicon

50.23.28

SPECIFICATIONS

[curved](#), [serrated longitudinal](#)