

## KOCHER-LANGENBECK RETRACTOR

40X11mm

21.5cm (8.6")

**SKU:** 064-130-40

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                    |
|-------------------|------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>               |
| <b>TYPE</b>       | <a href="#">Retractors</a>         |
| <b>Aesculap</b>   | BT359                              |
| <b>KLS Martin</b> | 15-317-05                          |
| <b>AUTHORS</b>    | <a href="#">KOCHER, LANGENBECK</a> |
| <b>Medicon</b>    | 23.25.06                           |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                    |
|-------------------|------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>               |
| <b>TYPE</b>       | <a href="#">Retractors</a>         |
| <b>Aesculap</b>   | BT359                              |
| <b>KLS Martin</b> | 15-317-05                          |
| <b>AUTHORS</b>    | <a href="#">KOCHER, LANGENBECK</a> |

**Medicon**

23.25.06