

## KOCHER-NIPPON HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

18.5cm (7.2")



**SKU:** 040-406-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-382-18
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">NIPPON</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>KLS Martin</b>	13-382-18
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">NIPPON</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>