

KOCHER-NIPPON HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

14.5cm (5.8")



SKU: 040-406-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-14
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-382-14
AUTHORS	KOCHER , NIPPON
SPECIFICATIONS	serrated , straight , toothed