

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

curved; serrated; 1X2 teeth

20.0cm (8")



SKU: 040-351-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH647
KLS Martin	13-415-20
AUTHORS	KOCHER , OCHSNER
Medicon	16.13.20
SPECIFICATIONS	curved , serrated , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH647
KLS Martin	13-415-20

AUTHORS

[KOCHER](#), [OCHSNER](#)

Medicon

16.13.20

SPECIFICATIONS

[curved](#), [serrated](#), [toothed](#)