

## KOCHER-OCHSNER HAEMOSTATIC FORCEPS

curved; serrated; 1X2 teeth

30.0cm (12")



**SKU:** 040-351-30

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH655
<b>KLS Martin</b>	13-415-30
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">OCHSNER</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH655
<b>KLS Martin</b>	13-415-30
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">OCHSNER</a>

## SPECIFICATIONS

[curved](#), [serrated](#), [toothed](#)