

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

16.0cm (6.4")



SKU: 040-350-16

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH642
KLS Martin	13-414-16
AUTHORS	KOCHER , OCHSNER
Medicon	16.12.16
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH642
KLS Martin	13-414-16

AUTHORS

[KOCHER](#), [OCHSNER](#)

Medicon

16.12.16

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)