

## KOCHER-OCHSNER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

14.0cm (5.6")



**SKU:** 040-350-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH610
<b>KLS Martin</b>	13-414-14
<b>AUTHORS</b>	<a href="#">KOCHER</a> , <a href="#">OCHSNER</a>
<b>Medicon</b>	16.12.14
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH610
<b>KLS Martin</b>	13-414-14

**AUTHORS**

[KOCHER](#), [OCHSNER](#)

**Medicon**

16.12.14

**SPECIFICATIONS**

[serrated](#), [straight](#), [toothed](#)