

KOCHER-OCHSNER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

14.0cm (5.6")



SKU: 040-350-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH610
KLS Martin	13-414-14
AUTHORS	KOCHER , OCHSNER
Medicon	16.12.14
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH610
KLS Martin	13-414-14

AUTHORS

[KOCHER](#), [OCHSNER](#)

Medicon

16.12.14

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)