

KOCHER RETRACTOR

3 teeth; semi-sharp; 15X15mm

22.0cm (8.8")


SKU: 064-013-03

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT263
KLS Martin	15-242-03
AUTHORS	KOCHER
Medicon	22.42.03
SPECIFICATIONS	semi-sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT263
KLS Martin	15-242-03

AUTHORS

[KOCHER](#)

Medicon

22.42.03

SPECIFICATIONS

[semi-sharp](#), [toothed](#)