

KÖNIG LIGATURE CONDUCTOR

curved; 5mm

19.5cm (7.8")

SKU: 056-508-05

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Conductors |
| Aesculap | BM901 |
| KLS Martin | 20-875-05 |
| AUTHORS | KOENIG |
| Medicon | 13.32.05 |
| SPECIFICATIONS | curved |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Ligature Conductors |
| Aesculap | BM901 |
| KLS Martin | 20-875-05 |

AUTHORS

[KOENIG](#)

Medicon

13.32.05

SPECIFICATIONS

[curved](#)