

KÖNIG LIGATURE CONDUCTOR

curved; 7mm

19.5cm (7.8")


SKU: 056-508-07

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM902
KLS Martin	20-875-07
AUTHORS	KOENIG
Medicon	13.32.07
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM902
KLS Martin	20-875-07

AUTHORS

[KOENIG](#)

Medicon

13.32.07

SPECIFICATIONS

[curved](#)