

## KÖNIG LIGATURE CONDUCTOR

curved; 7mm

19.5cm (7.8")


**SKU:** 056-508-07

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM902
<b>KLS Martin</b>	20-875-07
<b>AUTHORS</b>	<a href="#">KOENIG</a>
<b>Medicon</b>	13.32.07
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM902
<b>KLS Martin</b>	20-875-07

**AUTHORS**

[KOENIG](#)

**Medicon**

13.32.07

**SPECIFICATIONS**

[curved](#)