

KÖNIG LIGATURE CONDUCTOR

curved; 3mm

19.5cm (7.8")


SKU: 056-508-03

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM900
KLS Martin	20-875-03
AUTHORS	KOENIG
Medicon	13.32.03
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM900
KLS Martin	20-875-03

AUTHORS

[KOENIG](#)

Medicon

13.32.03

SPECIFICATIONS

[curved](#)