

## KRAUSE EXCHANGEABLE TIP

fig. 2; 3.5mm; through-cutting

for BRÜNINGS interchangeable tubes

**SKU:** 131-489-02

**Categories:** [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP223
<b>KLS Martin</b>	39-628-02
<b>AUTHORS</b>	<a href="#">KRAUSE</a>
<b>Medicon</b>	73.43.35



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP223
<b>KLS Martin</b>	39-628-02
<b>AUTHORS</b>	<a href="#">KRAUSE</a>

**Medicon**

73.43.35