

## KROENLEIN HAEMOSTATIC FORCEPS

straight; serrated

13.0cm (5.2")



**SKU:** 040-369-13

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Artery Forceps](#)

**Aesculap** BH300

**KLS Martin** 13-182-13

**AUTHORS** [KROENLEIN](#)

**Medicon** 16.06.43

**SPECIFICATIONS** [serrated](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Artery Forceps](#)

**Aesculap** BH300

**KLS Martin** 13-182-13

**AUTHORS**

[KROENLEIN](#)

**Medicon**

16.06.43

**SPECIFICATIONS**

[serrated](#), [straight](#)