

KRONECKER LIGATURE NEEDLE

right; with round handle

15.0cm (6")

SKU: 056-470-01

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Needles
Aesculap	BM801
KLS Martin	20-835-13
AUTHORS	KRONECKER
Medicon	13.11.14



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Needles
Aesculap	BM801
KLS Martin	20-835-13
AUTHORS	KRONECKER

Medicon

13.11.14