

## KRONECKER LIGATURE NEEDLE

right; with round handle

15.0cm (6")

**SKU:** 056-470-01

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Needles</a>
<b>Aesculap</b>	BM801
<b>KLS Martin</b>	20-835-13
<b>AUTHORS</b>	<a href="#">KRONECKER</a>
<b>Medicon</b>	13.11.14



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Needles</a>
<b>Aesculap</b>	BM801
<b>KLS Martin</b>	20-835-13
<b>AUTHORS</b>	<a href="#">KRONECKER</a>

**Medicon**

13.11.14