

## L-HOOK ELECTRODE; ADAPTABLE

for pistol handle 620-010-00

Ø 5.0mm; 330mm



**SKU:** 622-050-33

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes, monopolar](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [HF-Electrodes, monopolar](#)