

## LANE TOWEL FORCEPS

straight; toothed

20.0cm (8")




**SKU:** 030-040-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF452
<b>KLS Martin</b>	13-935-20
<b>AUTHORS</b>	<a href="#">LANE</a>
<b>Medicon</b>	17.59.20
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF452
<b>KLS Martin</b>	13-935-20

## AUTHORS

[LANE](#)

## Medicon

17.59.20

## SPECIFICATIONS

[straight](#), [toothed](#)