

LANE TOWEL FORCEPS

straight; toothed

15.0cm (6")



SKU: 030-040-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Towel Forceps](#)

Aesculap BF451

AUTHORS [LANE](#)

Medicon 17.59.15

SPECIFICATIONS [straight](#), [toothed](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Towel Forceps](#)

Aesculap BF451

AUTHORS [LANE](#)

Medicon 17.59.15

SPECIFICATIONS

[straight](#), [toothed](#)