

LANGE EXCHANGEABLE TIP

fig. 1; oval

for BRÜNINGS interchangeable tubes



SKU: 131-488-00

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Universal handles and tips |
| Aesculap | OP220 |
| KLS Martin | 39-624-01 |
| AUTHORS | LANGE |
| Medicon | 73.43.15 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Universal handles and tips |
| Aesculap | OP220 |
| KLS Martin | 39-624-01 |
| AUTHORS | LANGE |

Medicon

73.43.15