

LANGE-HOHMANN BONE LEVER

22mm; wide tip with groove

25.0cm (10")


SKU: 190-574-22

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK160
KLS Martin	23-549-22
AUTHORS	HOHMANN , LANGE



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK160
KLS Martin	23-549-22
AUTHORS	HOHMANN , LANGE