

## LAPAROSCOPE HD+

Ø 10.0mm; 45°; 330 mm

autoclavable

**SKU:** 601-010-40HD

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)



45°

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Laparoscopy](#)