

## LAPAROSCOPE HD+

Ø 10.0mm; 30°; 425 mm

autoclavable



**SKU:** 601-012-30HD

**Categories:** [Endoscopy](#), [Equipment](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Telescopes</a>  |
| <b>SCOPE</b> | <a href="#">Laparoscopy</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Telescopes</a>  |
| <b>SCOPE</b> | <a href="#">Laparoscopy</a> |