

LAPAROSCOPE HD+

Ø 10.0mm; 30°; 425 mm

autoclavable



SKU: 601-012-30HD

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Laparoscopy |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------|
| UNIT | pcs. |
| RANGE | Telescopes |
| SCOPE | Laparoscopy |